

## Product datasheet for **SR317222**

### ZNF619 Human siRNA Oligo Duplex (Locus ID 285267)

#### Product data:

|                                |                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>           | siRNA Oligo Duplexes                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity:</b>                 | HPLC purified                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Quality Control:</b>        | Tested by ESI-MS                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Sequences:</b>              | Available with shipment                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Stability:</b>              | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                                                                                                                                                                                        |
| <b># of transfections:</b>     | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                                                                                                                                                                                      |
| <b>Note:</b>                   | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                 | <u><a href="#">NM_001145082</a></u> , <u><a href="#">NM_001145083</a></u> , <u><a href="#">NM_001145093</a></u> , <u><a href="#">NM_001145094</a></u> , <u><a href="#">NM_173656</a></u> , <u><a href="#">NR_026746</a></u> ,<br><u><a href="#">NM_001363277</a></u>                                                                                                                                                        |
| <b>UniProt ID:</b>             | <u><a href="#">Q8N2I2</a></u>                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Components:</b>             | ZNF619 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 285267)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml                                                                                                                                                                     |
| <b>Summary:</b>                | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                                                                                                                                                               |
| <b>Performance Guaranteed:</b> | OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency. |

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).

