

Product datasheet for **SR317593**

KCNT2 Human siRNA Oligo Duplex (Locus ID 343450)

Product data:

| | |
|---------------------|--|
| Product Type: | siRNA Oligo Duplexes |
| Purity: | HPLC purified |
| Quality Control: | Tested by ESI-MS |
| Sequences: | Available with shipment |
| Stability: | One year from date of shipment when stored at -20°C. |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM). |
| Note: | Single siRNA duplex (10nmol) can be ordered. |
| RefSeq: | NM_001287819 , NM_001287820 , NM_198503 , NR_146057 , NR_146058 |
| UniProt ID: | Q6UVM3 |
| Synonyms: | EIEE57; KCa4.2; SLICK; SLO2.1 |
| Components: | KCNT2 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 343450) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml |
| Summary: | Outward rectifying potassium channel. Produces rapidly activating outward rectifier K(+) currents. Activated by high intracellular sodium and chloride levels (PubMed:14684870, PubMed:16687497, PubMed:29069600). Channel activity is inhibited by ATP and by inhalation anesthetics, such as isoflurane (PubMed:16687497) (By similarity). Inhibited upon stimulation of G-protein coupled receptors, such as CHRM1 and GRM1 (PubMed:16687497). [UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).