

## Product datasheet for **SR317938**

### VSTM5 Human siRNA Oligo Duplex (Locus ID 387804)

#### Product data:

|                     |  |
|---------------------|--|
| Product Type:       | siRNA Oligo Duplexes   |
| Purity:             | HPLC purified  |
| Quality Control:    | Tested by ESI-MS   |
| Sequences:          | Available with shipment  |
| Stability:          | One year from date of shipment when stored at -20°C.   |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| Note:               | Single siRNA duplex (10nmol) can be ordered.   |
| RefSeq:             | <a href="#">NM_001144871</a>   |
| UniProt ID:         | <a href="#">A8MXK1</a>   |
| Synonyms:           | C11orf90   |
| Components:         | VSTM5 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 387804)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml   |
| Summary:            | Cell adhesion-like membrane protein of the central nervous system (CNS) which modulates both the position and complexity of central neurons by altering their membrane morphology and dynamics. Involved in the formation of neuronal dendrites and protrusions including dendritic filopodia. In synaptogenesis, regulates synapse formation by altering dendritic spine morphology and actin distribution. Promotes formation of unstable neuronal spines such as thin and branched types. Regulates neuronal morphogenesis and migration during cortical development in the brain.[UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).