

Product datasheet for **SR318617**

TARP Human siRNA Oligo Duplex (Locus ID 445347)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001003799 , NM_001003806
Synonyms:	CD3G; TCRG; TCRGC1; TCRGC2; TCRGV
Components:	TARP (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 445347) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	In some non-lymphoid tissues, the unrearranged T cell receptor gamma (TRG@) locus is expressed. The resulting transcript contains a subset of the TRG@ gene segments and is shorter than TRG@ transcripts expressed in lymphoid tissues. This RefSeq record represents the unrearranged TRG@ locus transcript; the complete TRG@ locus is represented by the genomic RefSeq NG_001336. The transcript represented by this RefSeq has two open reading frames (ORFs) that encode different proteins. The downstream ORF is in the same frame as TRG@ and its protein product is similar to TRG@ proteins. The upstream ORF uses a different reading frame and encodes a novel protein. [provided by RefSeq, Jul 2008]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).