

## Product datasheet for **SR322601**

### **MBD4 Human siRNA Oligo Duplex (Locus ID 8930)**

#### **Product data:**

|                            |   |
|----------------------------|---|
| <b>Product Type:</b>       | siRNA Oligo Duplexes  |
| <b>Purity:</b>             | HPLC purified   |
| <b>Quality Control:</b>    | Tested by ESI-MS  |
| <b>Sequences:</b>          | Available with shipment   |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.  |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).  |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.  |
| <b>RefSeq:</b>             | <a href="#">NM_001276270</a> , <a href="#">NM_001276271</a> , <a href="#">NM_001276272</a> , <a href="#">NM_001276273</a> , <a href="#">NM_003925</a>   |
| <b>UniProt ID:</b>         | <a href="#">O95243</a>  |
| <b>Synonyms:</b>           | MED1  |
| <b>Components:</b>         | MBD4 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 8930)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml   |
| <b>Summary:</b>            | The protein encoded by this gene is a member of a family of nuclear proteins related by the presence of a methyl-CpG binding domain (MBD). These proteins are capable of binding specifically to methylated DNA, and some members can also repress transcription from methylated gene promoters. This protein contains an MBD domain at the N-terminus that functions both in binding to methylated DNA and in protein interactions and a C-terminal mismatch-specific glycosylase domain that is involved in DNA repair. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2013] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).