

Product datasheet for **SR322888**

RNF144 (RNF144A) Human siRNA Oligo Duplex (Locus ID 9781)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_014746 , NM_001349181 , NM_001349182 , NM_001349183 , NM_001349184 , NM_001349185 , NM_001349186 , NR_146069
UniProt ID:	P50876
Synonyms:	hUIP4; RNF144; UBCE7IP4
Components:	RNF144A (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 9781) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of a family of RING finger domain-containing E3 ubiquitin ligases that also includes parkin and parc. The expression of this gene is induced by DNA damage. The encoded protein interacts with the cytoplasmic DNA-dependent protein kinase, catalytic subunit (DNA-PKcs) and promotes its degradation through ubiquitination. The orthologous mouse protein has been shown to interact with a ubiquitin-conjugating enzyme involved in embryonic development. [provided by RefSeq, Mar 2017]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).