

Product datasheet for **SR400627**

Ufm1 Mouse siRNA Oligo Duplex (Locus ID 67890)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_026435
UniProt ID:	P61961
Synonyms:	1810045K17Rik; AI132708; AI463323; ENSMUSG00000074598; Gm10726
Components:	Ufm1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 67890) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Ubiquitin-like modifier which can be covalently attached via an isopeptide bond to substrate proteins as a monomer or a lysine-linked polymer (PubMed:21494687). The so-called ufmylation, requires the UFM1-activating E1 enzyme UBA5, the UFM1-conjugating E2 enzyme UFC1, and the UFM1-ligase E3 enzyme UFL1. This post-translational modification on lysine residues of proteins may play a crucial role in a number of cellular processes. TRIP4 ufmylation may for instance play a role in nuclear receptors-mediated transcription (By similarity). Other substrates may include DDRGK1 with which it may play a role in the cellular response to endoplasmic reticulum stress.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).