

Product datasheet for **SR400770**

Scgb3a2 Mouse siRNA Oligo Duplex (Locus ID 117158)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001289643 , NM_001289644 , NM_054038
UniProt ID:	Q920H1
Synonyms:	LuLeu1; Pnsp1; UGRP1; Utgrp1
Components:	Scgb3a2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 117158) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Secreted cytokine-like protein (By similarity). Binds to the scavenger receptor MARCO (By similarity). Can also bind to pathogens including the Gram-positive bacterium <i>L.monocytogenes</i> , the Gram-negative bacterium <i>P.aeruginosa</i> , and yeast (By similarity). Strongly inhibits phospholipase A2 (PLA2G1B) activity (PubMed:24213919). Seems to have anti-inflammatory effects in respiratory epithelium (PubMed:16456148, PubMed:25242865). Also has anti-fibrotic activity in lung (PubMed:24213919, PubMed:26559674). May play a role in fetal lung development and maturation (PubMed:18535256). Promotes branching morphogenesis during early stages of lung development (PubMed:18535256). In the pituitary, may inhibit production of follicle-stimulating hormone (FSH) and luteinizing hormone (LH) (PubMed:24514953).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).