

## Product datasheet for **SR401022**

### Ap4s1 Mouse siRNA Oligo Duplex (Locus ID 11782)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_021710</a>
UniProt ID:	<a href="#">Q9WVL1</a>
Synonyms:	AI314282
Components:	Ap4s1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 11782) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes the sigma subunit of the adaptor-related protein complex 4 which mediates intracellular membrane trafficking along the endocytic and secretory transport pathways. This complex contains four subunits, beta, epsilon, mu, and sigma, and belongs to a family of five adapter protein complexes, including three clathrin-associated complexes and two non clathrin-associated complexes, that localize to different intracellular compartments and mediate membrane vesicle trafficking using distinct pathways. In humans, loss-of-function mutations in this gene have been linked to specific adapter complex 4 deficiency disorders including hereditary spastic paraplegia. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jul 2016]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).