

Product datasheet for **SR402172**

Calca Mouse siRNA Oligo Duplex (Locus ID 12310)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001033954 , NM_001289444 , NM_007587
UniProt ID:	Q99JA0
Synonyms:	C; CA; Calc; Calc1; Cg; Cgrp; CGRP-1; CGRP1; Ct; Ctn
Components:	Calca (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12310) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes the peptide hormones calcitonin, calcitonin gene-related peptide (CGRP) and katalcalcin. Alternative splicing of the mRNA results in multiple variants that encode either calcitonin or CGRP preproteins. Post-translational processing of the calcitonin and CGRP propeptides results in either calcitonin and katalcalcin, or CGRP, respectively. Calcitonin and katalcalcin modulate calcium levels in the blood stream. CGRP can function as a vasodilator and play a role in the transmission of pain. The human homolog of CGRP was found to have antimicrobial activity. [provided by RefSeq, Mar 2015]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).