

## Product datasheet for **SR402685**

### Atp5g3 Mouse siRNA Oligo Duplex (Locus ID 228033)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001301721</a> , <a href="#">NM_001301722</a> , <a href="#">NM_175015</a>
UniProt ID:	<a href="#">P56384</a>
Synonyms:	6030447M23; Atp5mc3
Components:	Atp5g3 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 228033) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	The protein encoded by this gene is a subunit of mitochondrial membrane ATP synthase, the enzyme that catalyzes ATP synthesis during oxidative phosphorylation. This gene encodes subunit 9, which is present in multiple copies in the transmembrane part of the ATP synthase complex. Phenotype and gene expression profiles suggest correlations between this gene and alcoholism- and obesity-related phenotypes. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Sep 2014]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).