

Product datasheet for **SR402769**

Csn1s2b Mouse siRNA Oligo Duplex (Locus ID 12992)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001301334 , NM_001301335 , NM_009973
UniProt ID:	P02664
Synonyms:	AW987150; Cs; Csnd; Csne
Components:	Csn1s2b (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12992) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene is a member of the alpha-s2-like casein gene family, and this gene product is a calcium-sensitive casein. Members of this gene family are organized as a gene cluster that is conserved in its order, but with greater conservation amongst orthologs than paralogs. The protein encoded by this gene interacts with other casein proteins to form a micelle structure, and is a major source of protein in milk. This structure is important for the transport of calcium, phosphate, and protein. Alternative splicing results in multiple transcript variants encoding different protein isoforms. [provided by RefSeq, Aug 2014]


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).