

Product datasheet for **SR403106**

Cebpg Mouse siRNA Oligo Duplex (Locus ID 12611)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_009884
UniProt ID:	P53568
Synonyms:	C/EBPgamma; C/EBP[g]; GPE1-BP; Gpe1bp; Ig/EBP
Components:	Cebpg (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12611) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Transcription factor that binds to the promoter and the enhancer regions of target genes (PubMed:21602272). Binds to the promoter and the enhancer of the immunoglobulin heavy chain (PubMed:2121606). Binds to GPE1, a cis-acting element in the G-CSF gene promoter (PubMed:1709121). Binds to the enhancer element PRE-I (positive regulatory element-I) of the IL-4 gene (By similarity). Binds to the promoter and the enhancer of the alpha-1-fetoprotein gene (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).