

## Product datasheet for **SR403478**

### Trappc6a Mouse siRNA Oligo Duplex (Locus ID 67091)

#### Product data:

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | siRNA Oligo Duplexes                                                                                                                                                                                                                                                                                                                                                                                                        |
| Purity:                 | HPLC purified                                                                                                                                                                                                                                                                                                                                                                                                               |
| Quality Control:        | Tested by ESI-MS                                                                                                                                                                                                                                                                                                                                                                                                            |
| Sequences:              | Available with shipment                                                                                                                                                                                                                                                                                                                                                                                                     |
| Stability:              | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                                                                                                                                                                                        |
| # of transfections:     | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                                                                                                                                                                                      |
| Note:                   | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                                                                                                                                                                                                |
| RefSeq:                 | <a href="#">NM_001331181</a> , <a href="#">NM_025960</a> , <a href="#">NR_138582</a> , <a href="#">NR_138583</a> , <a href="#">NR_138584</a>                                                                                                                                                                                                                                                                                |
| UniProt ID:             | <a href="#">Q78XR0</a>                                                                                                                                                                                                                                                                                                                                                                                                      |
| Synonyms:               | 1810073E21Rik; 4930519D19Rik; AI480686; mhyp; TRS33                                                                                                                                                                                                                                                                                                                                                                         |
| Components:             | Trappc6a (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 67091)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml                                                                                                                                                                    |
| Summary:                | May play a role in vesicular transport during the biogenesis of melanosomes.<br>[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                                                                                                                             |
| Performance Guaranteed: | OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency. |

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)