

## Product datasheet for **SR404453**

### IL1f8 Mouse siRNA Oligo Duplex (Locus ID 69677)

#### Product data:

|                     |  |
|---------------------|--|
| Product Type:       | siRNA Oligo Duplexes   |
| Purity:             | HPLC purified  |
| Quality Control:    | Tested by ESI-MS   |
| Sequences:          | Available with shipment  |
| Stability:          | One year from date of shipment when stored at -20°C.   |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| Note:               | Single siRNA duplex (10nmol) can be ordered.   |
| RefSeq:             | <a href="#">NM_027163</a>  |
| UniProt ID:         | <a href="#">Q9D6Z6</a>   |
| Synonyms:           | 2310043N20Rik; IL36b   |
| Components:         | IL1f8 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 69677)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| Summary:            | Cytokine that binds to and signals through the IL1RL2/IL-36R receptor which in turn activates NF-kappa-B and MAPK signaling pathways in target cells linked to a pro-inflammatory response. Part of the IL-36 signaling system that is thought to be present in epithelial barriers and to take part in local inflammatory response; similar to the IL-1 system with which it shares the coreceptor IL1RAP. Stimulates production of interleukin-6 and interleukin-8 in synovial fibroblasts, articular chondrocytes and mature adipocytes. Induces expression of a number of antimicrobial peptides including beta-defensin 4 and beta-defensin 103 as well as a number of matrix metalloproteases (By similarity). Seems to be involved in skin inflammatory response by acting on keratinocytes, dendritic cells and indirectly on T-cells to drive tissue infiltration, cell maturation and cell proliferation. Induces the production of proinflammatory cytokines in bone marrow-derived dendritic cells (BMDCs), including IL-12, IL-1 beta, IL-6, TNF-alpha and IL-23, and activates p38 MAPK phosphorylation in BMDCs. Involved in dendritic cell maturation by stimulating the surface expression of CD80, CD86 and MHC class II. Induces the production of IFN-gamma, IL-4 and IL-17 by T-helper 1 (Th1) cells, cultured CD4(+) T-cells and splenocytes.[UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).