

## Product datasheet for **SR405147**

### **Psmb2 Mouse siRNA Oligo Duplex (Locus ID 26445)**

#### **Product data:**

|                            |  |
|----------------------------|--|
| <b>Product Type:</b>       | siRNA Oligo Duplexes   |
| <b>Purity:</b>             | HPLC purified  |
| <b>Quality Control:</b>    | Tested by ESI-MS   |
| <b>Sequences:</b>          | Available with shipment  |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.   |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.   |
| <b>RefSeq:</b>             | <a href="#">NM_011970</a>  |
| <b>UniProt ID:</b>         | <a href="#">Q9R1P3</a>   |
| <b>Synonyms:</b>           | AU045357; AW108089; C7-I; D4Wsu33e; HC7-I  |
| <b>Components:</b>         | Psmb2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 26445)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>            | Component of the 20S core proteasome complex involved in the proteolytic degradation of most intracellular proteins. This complex plays numerous essential roles within the cell by associating with different regulatory particles. Associated with two 19S regulatory particles, forms the 26S proteasome and thus participates in the ATP-dependent degradation of ubiquitinated proteins. The 26S proteasome plays a key role in the maintenance of protein homeostasis by removing misfolded or damaged proteins that could impair cellular functions, and by removing proteins whose functions are no longer required. Associated with the PA200 or PA28, the 20S proteasome mediates ubiquitin-independent protein degradation. This type of proteolysis is required in several pathways including spermatogenesis (20S-PA200 complex) or generation of a subset of MHC class I-presented antigenic peptides (20S-PA28 complex).[UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).