

## Product datasheet for **SR405175**

### Vgll2 Mouse siRNA Oligo Duplex (Locus ID 215031)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001300957</a> , <a href="#">NM_153786</a>
UniProt ID:	<a href="#">Q8BGW8</a>
Synonyms:	C130057C21Rik; Vgl; vgl-2; Vi; VIT; VITO-1; Vito1
Components:	Vgll2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 215031) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene is a member of the Vestigial-like (Vgl) gene family and is upregulated during muscle differentiation. The product of this gene interacts with and modifies the DNA-binding properties of the transcription factor, TEF-1, and is important for muscle tissue development. Reduced expression of this gene leads to a reduction in the terminal differentiation of muscle cells. Alternate splicing results in multiple protein isoforms. [provided by RefSeq, Jul 2014]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).