

Product datasheet for **SR406267**

Gpihbp1 Mouse siRNA Oligo Duplex (Locus ID 68453)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001347039 , NM_026730 , NM_001361294
UniProt ID:	Q9D1N2
Synonyms:	1110002J19Rik; GPI-HBP1
Components:	Gpihbp1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 68453) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Mediates the transport of lipoprotein lipase LPL from the basolateral to the apical surface of endothelial cells in capillaries (PubMed:20620994). Anchors LPL on the surface of endothelial cells in the lumen of blood capillaries (PubMed:20620994, PubMed:24726386, PubMed:27811232). Thereby, plays an important role in lipolytic processing of chylomicrons by LPL, triglyceride metabolism and lipid homeostasis (PubMed:17403372). Binds chylomicrons and phospholipid particles that contain APOA5 (PubMed:18340083). Binds high-density lipoprotein (HDL) and plays a role in the uptake of lipids from HDL (PubMed:12496272).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).