

Product datasheet for **SR406276**

Eif4h Mouse siRNA Oligo Duplex (Locus ID 22384)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_033561
UniProt ID:	Q9WUK2
Synonyms:	AU018978; D5Ertd355; D5Ertd355e; E430026L18Rik; Ef4h; eIF-4H; Eif4; Wbsc; Wbscr1; Wsc; Wscr1
Components:	Eif4h (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 22384) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes eukaryotic translation initiation factor 4H (eIF4H) that plays a critical role in the process of protein synthesis. The encoded protein is an RNA-binding protein that, in concert with other translation initiation factors, helps unwind the 5' cap-proximal region of mRNA to prepare it for ribosomal attachment. Mice lacking the encoded protein displayed growth retardation with a significant reduction of body weight, a smaller brain volume and altered brain morphology. Behaviorally, mice lacking the encoded protein exhibit severe impairments of fear-related associative learning and memory formation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2015]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).