

Product datasheet for **SR406759**

Pbx1 Mouse siRNA Oligo Duplex (Locus ID 18514)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001291508 , NM_001291509 , NM_008783 , NM_183355
UniProt ID:	P41778
Synonyms:	2310056B04Rik; D230003C07Rik; Pbx; Pbx-; Pbx-1
Components:	Pbx1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 18514) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a homeobox protein that belongs to the three-amino-acid loop extension/Pre-B cell leukemia transcription factor (TALE/PBX) family of proteins. The encoded protein is involved in several biological processes during embryogenesis including steroidogenesis, sexual development and the maintenance of hematopoietic stem cells. This protein functions in the development of several organ systems and plays a role in skeletal patterning and programming. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).