

Product datasheet for **SR407293**

Casp14 Mouse siRNA Oligo Duplex (Locus ID 12365)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_009809
UniProt ID:	O89094
Synonyms:	MICE; mini-ICE
Components:	Casp14 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12365) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Non-apoptotic caspase which is involved in epidermal differentiation. Seems to play a role in keratinocyte differentiation and is required for cornification (PubMed:18156206). Regulates maturation of the epidermis by proteolytically processing filaggrin (PubMed:21654840). In vitro is equally active on the synthetic caspase substrates WEHD-ACF and IETD-AFC. Involved in processing of prosaposin in the epidermis (PubMed:24872419). May be involved in retinal pigment epithelium cell barrier function (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).