

Product datasheet for **SR407349**

Emd Mouse siRNA Oligo Duplex (Locus ID 13726)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_007927
UniProt ID:	O08579
Synonyms:	AW550900; Sta
Components:	Emd (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 13726) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene is the founding member of the LEM (Lap2beta, Emerin, and Man1) domain gene family. It encodes an integral membrane protein of the inner nuclear membrane. This gene may be involved in transcriptional regulation and genomic organization. In humans, mutations in this gene cause a class of diseases called laminopathies, which includes Emery-Dreifuss muscular dystrophy (EDMD). [provided by RefSeq, Apr 2013]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).