

## Product datasheet for **SR407565**

### Gzmm Mouse siRNA Oligo Duplex (Locus ID 16904)

#### Product data:

|                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:       | siRNA Oligo Duplexes                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Purity:             | HPLC purified                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Quality Control:    | Tested by ESI-MS                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Sequences:          | Available with shipment                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Stability:          | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Note:               | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| RefSeq:             | <a href="#">NM_001302485</a> , <a href="#">NM_001302499</a> , <a href="#">NM_008504</a> , <a href="#">NR_126325</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Synonyms:           | Lmet; Lmet1; MMET-1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Components:         | Gzmm (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 16904)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml                                                                                                                                                                                                                                                                                                                 |
| Summary:            | The protein encoded by this gene is a member of a family of cytotoxic lymphocyte serine proteases called granzymes, which are expressed by cytotoxic T lymphocytes and natural killer cells. This protein belongs to a subfamily of granzymes that cleave after methionine residues. Natural killer cell development, homeostasis and cytotoxicity are normal in mice deficient for this gene, but they demonstrate increased susceptibility to murine cytomegalovirus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).