

Product datasheet for **SR407784**

Cela2a Mouse siRNA Oligo Duplex (Locus ID 13706)

Product data:

| | |
|---------------------|--|
| Product Type: | siRNA Oligo Duplexes |
| Purity: | HPLC purified |
| Quality Control: | Tested by ESI-MS |
| Sequences: | Available with shipment |
| Stability: | One year from date of shipment when stored at -20°C. |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM). |
| Note: | Single siRNA duplex (10nmol) can be ordered. |
| RefSeq: | NM_007919 |
| UniProt ID: | P05208 |
| Synonyms: | EI; Ela; Ela-2; Ela2; Ela2a |
| Components: | Cela2a (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 13706) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml |
| Summary: | This gene encodes a serine protease enzyme that hydrolyzes elastin. This gene is highly expressed in the pancreatic acinar cells where the encoded preproprotein undergoes processing including signal peptide cleavage to generate an inactive zymogen. The removal of N-terminal activation peptide from the zymogen by trypsin generates active elastase enzyme. This gene is also expressed in the mouse epidermis where it participates in pro-filaggrin processing. [provided by RefSeq, Jul 2016] |



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).