

Product datasheet for **SR408189**

Bcdin3d Mouse siRNA Oligo Duplex (Locus ID 75284)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_029236
UniProt ID:	Q91YP1
Synonyms:	4930556P03Rik; AV138748
Components:	Bcdin3d (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 75284) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	O-methyltransferase that specifically monomethylates 5'-monophosphate of cytoplasmic histidyl tRNA, acting as a capping enzyme. Less efficiently, also methylates the 5' monophosphate of pre-miRNAs, acting as a negative regulator of miRNA processing. The 5' monophosphate of pre-miRNAs is recognized by DICER1 and is required for pre-miRNAs processing: methylation at this position reduces the processing of pre-miRNAs by DICER1. Able to mediate methylation of pre-miR-145, as well as other pre-miRNAs. There is some controversy about the methylation of pre-miR-145, since the dimethylation first described as the specific enzymatic activity cannot be reproduced by a more recent work which observes a monomethylation of pre-miR-145 but two orders weaker than the methylation of cytosolic histidyl tRNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).