

Product datasheet for **SR408849**

Medag Mouse siRNA Oligo Duplex (Locus ID 70717)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<u>NM_027519</u> , <u>NM_001356995</u> , <u>NR_151595</u>
UniProt ID:	<u>Q14BA6</u>
Synonyms:	6330406I15Rik; Awms3; mAWMS3; MEDA-4; Meda4
Components:	Medag (Mouse) – 3 unique 27mer siRNA duplexes – 2 nmol each (Locus ID 70717) Included – SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex – 2 nmol Included – SR30005, RNase free siRNA Duplex Resuspension Buffer – 2 ml
Summary:	Involved in processes that promote adipocyte differentiation, lipid accumulation, and glucose uptake in mature adipocytes.[UniProtKB/Swiss-Prot Function]



Performance Guaranteed: OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).