

## Product datasheet for **SR409888**

### Hoxa11 Mouse siRNA Oligo Duplex (Locus ID 15396)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_010450</a>
UniProt ID:	<a href="#">P31311</a>
Synonyms:	Hox-1.; Hox-1.9; Hoxa-; Hoxa-11
Components:	Hoxa11 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 15396) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene is located in a cluster of developmentally and temporally regulated genes on chromosome 6 encoding proteins involved in pattern formation. These proteins contain a characteristic DNA-binding motif called a homeodomain and function in transcriptional regulation. There are four distinct clusters of related genes on chromosomes 2, 6, 11, and 15. The protein encoded by this gene is important in the development of the skeleton, limbs, and urogenital tract. Expression of this gene may be regulated by overlapping transcription from an adjacent locus on the opposite strand (GeneID: 15397). [provided by RefSeq, Mar 2013]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).