

Product datasheet for SR410262

Rdh11 Mouse siRNA Oligo Duplex (Locus ID 17252)

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<u>NM 021557, NM 001362269, NM 001362270</u>
UniProt ID:	Q9QYF1
Synonyms:	2610319N22Rik; Al428145; Arsdr1; AU045252; C85936; CGI-82; HCBP12; M42C60; Mdt1; Psdr1; SCALD
Components:	Rdh11 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 17252) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Retinol dehydrogenase with a clear preference for NADP (PubMed:12807874, PubMed:29567832). Displays high activity towards 9-cis, 11-cis and all-trans-retinol, and to a lesser extent on 13-cis-retinol (By similarity) (PubMed:12807874). Exhibits also reductive activity towards toxic lipid peroxidation products such as medium-chain aldehydes trans-2- nonenal, nonanal, and cis-6-nonenal (PubMed:12807874). Has no dehydrogenase activity towards steroid (PubMed:12807874). Seems to be required for homeostasis of retinol in liver and testis (PubMed:29567832).[UniProtKB/Swiss-Prot Function]



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2021 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Rdh11 Mouse siRNA Oligo Duplex (Locus ID 17252) – SR410262
Performance Guaranteed:	OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.
	For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2021 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US