

Product datasheet for **SR411344**

Septin3 Mouse siRNA Oligo Duplex (Locus ID 24050)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_011889 , NM_001358836 , NM_001358837 , NM_001368642 , NM_001368643 , NM_001368644 , NM_001368645 , NM_001368646 , NM_001368647 , NM_001368784 , NM_001368785
UniProt ID:	Q9Z1S5
Synonyms:	3110018K01Rik; AV154067; B530002E20Rik; Sep3
Components:	Sept3 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 24050) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Filament-forming cytoskeletal GTPase (By similarity). May play a role in cytokinesis (Potential). [UniProtKB/Swiss-Prot Function]
Performance Guaranteed:	OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)