

## Product datasheet for **SR412048**

### Ptgdr Mouse siRNA Oligo Duplex (Locus ID 19214)

#### Product data:

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | siRNA Oligo Duplexes                                                                                                                                                                                                                                                                                                                                                                                                        |
| Purity:                 | HPLC purified                                                                                                                                                                                                                                                                                                                                                                                                               |
| Quality Control:        | Tested by ESI-MS                                                                                                                                                                                                                                                                                                                                                                                                            |
| Sequences:              | Available with shipment                                                                                                                                                                                                                                                                                                                                                                                                     |
| Stability:              | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                                                                                                                                                                                        |
| # of transfections:     | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                                                                                                                                                                                      |
| Note:                   | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                                                                                                                                                                                                |
| RefSeq:                 | <a href="#">NM_008962</a>                                                                                                                                                                                                                                                                                                                                                                                                   |
| UniProt ID:             | <a href="#">P70263</a>                                                                                                                                                                                                                                                                                                                                                                                                      |
| Synonyms:               | DP; PGD                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Components:             | Ptgdr (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 19214)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml                                                                                                                                                                       |
| Summary:                | Receptor for prostaglandin D2 (PGD2). The activity of this receptor is mainly mediated by G(s) proteins that stimulate adenylate cyclase, resulting in an elevation of intracellular cAMP. A mobilization of calcium is also observed, but without formation of inositol 1,4,5-trisphosphate (By similarity).[UniProtKB/Swiss-Prot Function]                                                                                |
| Performance Guaranteed: | OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency. |

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)