

## Product datasheet for **SR412598**

### Glul Mouse siRNA Oligo Duplex (Locus ID 14645)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_008131</a>
UniProt ID:	<a href="#">P15105</a>
Synonyms:	Glns; GS
Components:	Glul (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 14645) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Glutamine synthetase that catalyzes the ATP-dependent conversion of glutamate and ammonia to glutamine (By similarity). Its role depends on tissue localization: in the brain, it regulates the levels of toxic ammonia and converts neurotoxic glutamate to harmless glutamine, whereas in the liver, it is one of the enzymes responsible for the removal of ammonia (PubMed:25870278). Essential for proliferation of fetal skin fibroblasts (By similarity). Independently of its glutamine synthetase activity, required for endothelial cell migration during vascular development (PubMed:30158707). Involved in angiogenesis by regulating membrane localization and activation of the GTPase RHOJ, possibly by promoting RHOJ palmitoylation (By similarity). May act as a palmitoyltransferase for RHOJ: able to autopalmitoylate and then transfer the palmitoyl group to RHOJ (By similarity). Plays a role in ribosomal 40S subunit biogenesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).