

Product datasheet for **SR412922**

Ptges2 Mouse siRNA Oligo Duplex (Locus ID 96979)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_133783
UniProt ID:	Q8BWM0
Synonyms:	0610038H10Rik; C79137; Gbf1; Mptes2; Pges2
Components:	Ptges2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 96979) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	The protein encoded by this gene is a Golgi membrane-associated prostaglandin E synthase candidate, which is capable of catalyzing the conversion of prostaglandin H2 to prostaglandin E2 in vitro. However, a study using mice deficient of this gene suggests that this enzyme does not contribute to prostaglandin E2 biosynthesis in vivo. This protein is synthesized as a Golgi membrane-bound protein, but its N-terminal hydrophobic region is cleaved off during protein maturation to produce the predominant soluble truncated form that still retains the enzyme activity. This soluble protein also has been shown to activate the transcription regulated by a gamma-interferon-activated transcription element (GATE), possibly via an interaction with CAAAT/enhancer-binding protein-beta. [provided by RefSeq, Oct 2009]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).