

Product datasheet for **SR412985**

Adipor2 Mouse siRNA Oligo Duplex (Locus ID 68465)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_197985 , NM_001355692
UniProt ID:	Q8BQS5
Synonyms:	1110001114Rik; ADCR2; AI115388; AW554121; D6Ucla1e; Paqr2
Components:	Adipor2 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 68465) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Receptor for ADIPOQ, an essential hormone secreted by adipocytes that regulates glucose and lipid metabolism (PubMed:17327425, PubMed:17068142, PubMed:17268472, PubMed:24742672). Required for normal body fat and glucose homeostasis (PubMed:17327425, PubMed:17068142, PubMed:17268472, PubMed:24742672). ADIPOQ-binding activates a signaling cascade that leads to increased PPARA activity, and ultimately to increased fatty acid oxidation and glucose uptake (PubMed:12802337, PubMed:17268472, PubMed:24742672). Has intermediate affinity for globular and full-length adiponectin (PubMed:12802337). Required for normal revascularization after chronic ischemia caused by severing of blood vessels (PubMed:24742672).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).