

## Product datasheet for **SR414155**

### Pramel7 Mouse siRNA Oligo Duplex (Locus ID 347712)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<u><a href="#">NM_178250</a></u>
UniProt ID:	<u><a href="#">Q810Y8</a></u>
Synonyms:	AA675346; C81180
Components:	Pramel7 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 347712) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Promotes maintenance and self-renewal of pluripotent embryonic stem cells (ESCs), downstream of LIF/STAT3 (PubMed:21425410). Maintains the pluripotency state of ESCs by repressing DNA methylation through the regulation of UHRF1 stability. Mediates the proteasomal degradation of UHRF1. Is required for the establishment of the blastocyst (PubMed:28604677).[UniProtKB/Swiss-Prot Function]

[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).