

## Product datasheet for **SR414632**

### Tph1 Mouse siRNA Oligo Duplex (Locus ID 21990)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001136084</a> , <a href="#">NM_001276372</a> , <a href="#">NM_009414</a>
UniProt ID:	<a href="#">P17532</a>
Synonyms:	Tph
Components:	Tph1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 21990) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the bipterin-dependent aromatic amino acid hydroxylase family. The encoded protein is one of two tryptophan hydroxylase enzymes that catalyze the first and rate limiting step in the biosynthesis of the hormone and neurotransmitter, serotonin. This gene is expressed in peripheral organs, while tryptophan hydroxylase 2 is expressed in neurons. The encoded protein is involved in the development of hypoxia-induced elevations in pulmonary pressures and pulmonary vascular remodeling, and has also been implicated as a regulator of immune tolerance. Disruption of this gene is associated with cardiac dysfunction. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2013]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).