

## Product datasheet for **SR414759**

### Cckbr Mouse siRNA Oligo Duplex (Locus ID 12426)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_007627</a>
UniProt ID:	<a href="#">P56481</a>
Synonyms:	CCK-BR; CCK2; CCK2-R; CCK2R; CCKR-; CCKR-2
Components:	Cckbr (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12426) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a multipass transmembrane receptor protein expressed in the central nervous system and gastrointestinal tract. Cholecystokinin and gastrin bind to the encoded protein to stimulate gastric acid secretion and mucosal growth in the gastrointestinal tract, and anxiety, pain sensation and memory in the brain. Mice lacking the encoded protein exhibit an increase in the basal gastric pH and gastrin levels in the bloodstream as well as mild hypocalcemia, secondary hyperparathyroidism and increased bone resorption. [provided by RefSeq, Apr 2015]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).