

Product datasheet for **SR415473**

Cpvl Mouse siRNA Oligo Duplex (Locus ID 71287)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001289713 , NM_001289714 , NM_027749
UniProt ID:	Q9D3S9
Synonyms:	4933436L16Rik; HVLP
Components:	Cpvl (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 71287) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the serine carboxypeptidase family of proteases that cleave amino acids from the C-terminus of a protein substrate. The human ortholog of this gene, where it was first characterized, was found to be upregulated during the maturation of monocytes to macrophages. The encoded protein may be involved in antigen processing, digestion of phagocytosed proteins in the lysosome and lamellipodium formation. Disruption of this gene in mice was found to cause embryonic lethality. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2015]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).