

## Product datasheet for **SR415545**

### **Fgb Mouse siRNA Oligo Duplex (Locus ID 110135)**

#### **Product data:**

<b>Product Type:</b>	siRNA Oligo Duplexes
<b>Purity:</b>	HPLC purified
<b>Quality Control:</b>	Tested by ESI-MS
<b>Sequences:</b>	Available with shipment
<b>Stability:</b>	One year from date of shipment when stored at -20°C.
<b># of transfections:</b>	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
<b>Note:</b>	Single siRNA duplex (10nmol) can be ordered.
<b>RefSeq:</b>	<a href="#">NM_181849</a>
<b>UniProt ID:</b>	<a href="#">Q8K0E8</a>
<b>Synonyms:</b>	2510049G14Rik
<b>Components:</b>	Fgb (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 110135) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
<b>Summary:</b>	This gene encodes the beta subunit of the coagulation factor fibrinogen, which is a component of the blood clot. The encoded preproprotein is proteolytically processed by thrombin to release an N-terminal fibrinopeptide during the conversion of fibrinogen to insoluble fibrin polymer. The encoded protein interacts with the amyloid beta peptide to form fibrin clots of abnormal structure, and may play an important role in Alzheimer's disease. This gene is located adjacent to the genes encoding fibrinogen alpha and gamma subunits on chromosome 3. [provided by RefSeq, Nov 2015]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).