

## Product datasheet for **SR415685**

### Dclre1c Mouse siRNA Oligo Duplex (Locus ID 227525)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001110214</a> , <a href="#">NM_001302674</a> , <a href="#">NM_001302684</a> , <a href="#">NM_146114</a> , <a href="#">NM_175683</a>
UniProt ID:	<a href="#">Q8K4J0</a>
Synonyms:	9930121L06Rik; A; AI661365; Art; Snm1l
Components:	Dclre1c (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 227525) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the SNM1 family of nucleases and is involved in V(D)J recombination and DNA repair. This protein has single-strand-specific 5'-3' exonuclease activity; it also exhibits endonuclease activity on 5' and 3' overhangs and hairpins. The protein also functions in the regulation of the cell cycle in response to DNA damage. Homozygous knockout mice for this gene exhibit severe combined immunodeficiency with sensitivity to ionizing radiation. Mutations in this gene in humans can cause Athabaskan-type severe combined immunodeficiency (SCIDA) and Omenn syndrome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2014]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).