

Product datasheet for **SR416137**

Bmpr1b Mouse siRNA Oligo Duplex (Locus ID 12167)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001277216 , NM_001277217 , NM_001277218 , NM_001277220 , NM_007560
UniProt ID:	P36898
Synonyms:	Acv; Acvrlk6; AI385617; AI; ALK-6; Alk6; AV355320; BMPR; BMPR-1B; BMPR-IB; CFK-43; CFK-43a; SKR6
Components:	Bmpr1b (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 12167) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a serine/threonine kinase that functions as a receptor for bone morphogenetic proteins (BMPs). The encoded protein is a type I receptor, and forms a complex of two type II and two type I receptors at the cell membrane. This complex signals downstream to activate SMAD transcriptional regulators. This signaling is important in skeletal and bone development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).