

## Product datasheet for **SR418517**

### Ano3 Mouse siRNA Oligo Duplex (Locus ID 228432)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001081556</a> , <a href="#">NM_001128103</a> , <a href="#">NM_001355244</a> , <a href="#">NM_001355245</a> , <a href="#">NM_001355246</a>
UniProt ID:	<a href="#">A2AHL1</a>
Synonyms:	A1838058; B230324K02Rik; Tmem16c
Components:	Ano3 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 228432) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Has calcium-dependent phospholipid scramblase activity; scrambles phosphatidylcholine and galactosylceramide (PubMed:23532839). Does not exhibit calcium-activated chloride channel (CaCC) activity (PubMed:23532839). Seems to act as potassium channel regulator and may inhibit pain signaling; can facilitate KCNT1/Slack channel activity by promoting its full single-channel conductance at very low sodium concentrations and by increasing its sodium sensitivity (PubMed:23872594).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).