

Product datasheet for **SR418927**

F13b Mouse siRNA Oligo Duplex (Locus ID 14060)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_031164
UniProt ID:	Q07968
Synonyms:	Cf-13b; Cf13b
Components:	F13b (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 14060) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes subunit B of the coagulation factor XIII that catalyzes the final step of the blood coagulation pathway. The encoded protein associates with subunit A to form a heterotetrameric protein that circulates in the plasma. During the blood coagulation process, thrombin-mediated proteolytic cleavage of the subunit A results in the dissociation of the encoded protein from the heterotetramer. Male mice lacking the encoded protein exhibit mild fibrosis together with hemosiderin deposits in the heart. [provided by RefSeq, Sep 2015]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).