

## Product datasheet for **SR419111**

### **E4f1 Mouse siRNA Oligo Duplex (Locus ID 13560)**

#### **Product data:**

|                                |   |
|--------------------------------|---|
| <b>Product Type:</b>           | siRNA Oligo Duplexes  |
| <b>Purity:</b>                 | HPLC purified   |
| <b>Quality Control:</b>        | Tested by ESI-MS  |
| <b>Sequences:</b>              | Available with shipment   |
| <b>Stability:</b>              | One year from date of shipment when stored at -20°C.  |
| <b># of transfections:</b>     | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).  |
| <b>Note:</b>                   | Single siRNA duplex (10nmol) can be ordered.  |
| <b>RefSeq:</b>                 | <a href="#">NM_001301784</a> , <a href="#">NM_007893</a> , <a href="#">NR_125990</a>  |
| <b>UniProt ID:</b>             | <a href="#">Q8CCE9</a>  |
| <b>Synonyms:</b>               | p120E4F; phi-AP3  |
| <b>Components:</b>             | E4f1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 13560)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>                | This gene encodes a member of the GLI-Kruppel zinc finger family. The encoded protein is likely to be multi-functional, with both adenovirus E1A-regulated transcription factor and ubiquitin E3 ligase activities, including roles in cell cycle regulation and the ubiquitination of p53. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2014]                                    |
| <b>Performance Guaranteed:</b> | OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency. |

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).


[View online »](#)