

## Product datasheet for **SR419440**

### Nckipsd Mouse siRNA Oligo Duplex (Locus ID 80987)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_030729</a>
UniProt ID:	<a href="#">Q9ESJ4</a>
Synonyms:	AF3P21; DIP1; ORF1; SPIN90; Wasbp; WASLBP; WISH
Components:	Nckipsd (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 80987) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Has an important role in stress fiber formation induced by active diaphanous protein homolog 1 (DRF1) (By similarity). Induces microspike formation, in vivo. In vitro, stimulates N-WASP-induced ARP2/3 complex activation in the absence of CDC42. May play an important role in the maintenance of sarcomere and/or in the assembly of myofibrils into sarcomeres. Implicated in regulation of actin polymerization and cell adhesion.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).