

## Product datasheet for **SR419795**

### Pla2g4a Mouse siRNA Oligo Duplex (Locus ID 18783)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_008869</a>
UniProt ID:	<a href="#">P47713</a>
Synonyms:	cP; cPL; cPLA2; cPLA2-alpha; cPLA2alpha; Pla; Pla2g4
Components:	Pla2g4a (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 18783) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	The protein encoded by this gene is a member of the phospholipase A2 group IV family. This enzyme hydrolyzes membrane phospholipids, thereby releasing the polyunsaturated fatty acid, arachidonic acid. Arachidonic acid is further metabolized into eicosanoids such as leukotrienes, thromboxanes and prostaglandins, that play important roles in regulating diverse biological processes such as inflammatory responses, membrane and actin dynamics, and tumorigenesis. A rise in intracellular calcium levels results in binding of calcium to the C2 domain of this protein, and triggers the translocation from the cytosol to intracellular membranes, including the Golgi apparatus. Disruption of this gene in mice led to decreased levels of eicosanoids and platelet-activating factor, decreased allergic symptoms, and impaired reproductive ability in females. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2015]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).