

## Product datasheet for **SR420355**

### Eml1 Mouse siRNA Oligo Duplex (Locus ID 68519)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001043335</a> , <a href="#">NM_001043336</a> , <a href="#">NM_001286346</a> , <a href="#">NM_001286347</a> , <a href="#">NM_001364193</a> , <a href="#">NM_001364194</a> , <a href="#">NM_001364195</a> , <a href="#">NM_001364196</a> , <a href="#">NM_001364197</a> , <a href="#">NM_001364198</a> , <a href="#">NM_001364199</a>
UniProt ID:	<a href="#">Q05BC3</a>
Synonyms:	1110008N23Rik; A930030P13Rik; AA171013; AI847476; AI853955; ELP79; EMAP; EMAPL; heco
Components:	Eml1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 68519) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Modulates the assembly and organization of the microtubule cytoskeleton, and probably plays a role in regulating the orientation of the mitotic spindle and the orientation of the plane of cell division. Required for normal proliferation of neuronal progenitor cells in the developing brain and for normal brain development. Does not affect neuron migration per se.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).