

## Product datasheet for **SR420496**

### **Cdhr5 Mouse siRNA Oligo Duplex (Locus ID 72040)**

#### **Product data:**

|                            |  |
|----------------------------|--|
| <b>Product Type:</b>       | siRNA Oligo Duplexes   |
| <b>Purity:</b>             | HPLC purified  |
| <b>Quality Control:</b>    | Tested by ESI-MS   |
| <b>Sequences:</b>          | Available with shipment  |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.   |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.   |
| <b>RefSeq:</b>             | <a href="#">NM_001114322</a> , <a href="#">NM_028069</a>   |
| <b>UniProt ID:</b>         | <a href="#">Q8VHF2</a>   |
| <b>Synonyms:</b>           | 1810074H01Rik; AI481143; Mucdhl; Mupcdh  |
| <b>Components:</b>         | Cdhr5 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 72040)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>            | Intermicrovillar adhesion molecule that forms, via its extracellular domain, calcium-dependent heterophilic complexes with CDHR2 on adjacent microvilli. Thereby, controls the packing of microvilli at the apical membrane of epithelial cells. Through its cytoplasmic domain, interacts with microvillus cytoplasmic proteins to form the intermicrovillar adhesion complex/IMAC. This complex plays a central role in microvilli and epithelial brush border differentiation.[UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).