

## Product datasheet for **SR421166**

### Hkdc1 Mouse siRNA Oligo Duplex (Locus ID 216019)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_145419</a>
UniProt ID:	<a href="#">Q91W97</a>
Synonyms:	BC016235
Components:	Hkdc1 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 216019) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Catalyzes the phosphorylation of hexose to hexose 6-phosphate, although at very low level compared to other hexokinases (By similarity). Has low glucose phosphorylating activity compared to other hexokinases (By similarity). Involved in glucose homeostasis and hepatic lipid accumulation (PubMed:30543855). Required to maintain whole-body glucose homeostasis during pregnancy; however additional evidences are required to confirm this role (PubMed:25648650, PubMed:27459389).[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).