

## Product datasheet for **SR421481**

### **Pcdh8 Mouse siRNA Oligo Duplex (Locus ID 18530)**

#### **Product data:**

|                            |   |
|----------------------------|---|
| <b>Product Type:</b>       | siRNA Oligo Duplexes  |
| <b>Purity:</b>             | HPLC purified   |
| <b>Quality Control:</b>    | Tested by ESI-MS  |
| <b>Sequences:</b>          | Available with shipment   |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.  |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).  |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.  |
| <b>RefSeq:</b>             | <a href="#">NM_001042726</a> , <a href="#">NM_021543</a>  |
| <b>UniProt ID:</b>         | <a href="#">Q7TSK3</a>  |
| <b>Synonyms:</b>           | 1700080P15Rik; P; Papc  |
| <b>Components:</b>         | Pcdh8 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 18530)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml   |
| <b>Summary:</b>            | This gene belongs to the protocadherin gene family, a subfamily of the cadherin superfamily. The gene encodes a type I transmembrane protein composed of an extracellular domain including 6 cadherin ectodomains, a single-pass transmembrane domain and a cytoplasmic tail. Unlike classical cadherins, which are generally encoded by 15-17 exons, this gene includes only 3 exons with the first large exon encoding the extracellular and transmembrane region. Although this gene product is capable of homophilic interaction, it appears to affect cell-cell adhesion indirectly by initiating signaling events that regulate classical cadherin-mediated adhesion. Based on studies on this protein and its orthologs, this protocadherin mainly functions in developing embryos and the central nervous system, but can also function as a tumor suppressor. Alternative splicing yielding isoforms with unique cytoplasmic tails has been reported. [provided by RefSeq, Sep 2009] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).