

Product datasheet for **SR421835**

Dsg4 Mouse siRNA Oligo Duplex (Locus ID 16769)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_181564
UniProt ID:	Q7TMD7
Synonyms:	CDHF13; I; lah
Components:	Dsg4 (Mouse) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 16769) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the cadherin family of proteins that forms an integral transmembrane component of desmosomes, the multiprotein complexes involved in cell adhesion, organization of the cytoskeleton, cell sorting and cell signaling. This gene is expressed in the suprabasal epidermis and hair follicle. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional protein. Certain mutations in this gene are responsible for the lanceolate hair phenotype in mice. This gene is located in a cluster of desmosomal cadherin genes on chromosome 18. [provided by RefSeq, Feb 2016]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).